

INTRODUCING
The **NEW** Omnigon
Support Garment

SUPPORT LEVEL

**Rating 3 -
MODERATE
SUPPORT**

SUITABLE FOR

- After surgery
- Wear during light exercise
- To help prevent or to support a small bulge

KomfortLite Support Belt

INSIDE POCKET

OUTSIDE POCKET

With handy pockets to help you grip the belt, hold it against your body, and pull the belt together.



OMNIGON
SUPPORT GARMENTS

Designed for comfort and ease-of-use

Soft breathable fabric for **comfort against your skin**

Available in black, with either a **17cm or 21.5cm width**

Designed with usability in mind. The inside hand pocket is used to help hold the belt in place, whilst the outside hand pocket helps bring the two ends around your waist.



Can be worn during **light exercise** – including swimming!

Strong hook and loop closure ensures belt is **secure and stays in place**

Rounded corners for flexible application, no matter what your body shape is. You can apply the belt at an angle, and/or may find it comfortable to wear the joining pads along your back.

Other benefits

- ✔ Disguise your stoma and pouch when worn underneath your shirt
- ✔ Smooth your body contour to give you the confidence during the day
- ✔ Provide pouch discretion
- ✔ Minimise noises and muffle sounds
- ✔ Draw in and help support the weight of bulges
- ✔ Feel more secure in your core

“

The comfort level for the KomfortLite support belt was beautiful, the fabric is very soft and felt smooth and delicate on the skin.

Wearing it in water was cool and allows you to have extra support and didn't become too heavy and dried very quickly.

A great product that will be very useful for so many ostomates. ”

JESS
@ownyourbag



Find the right size for you

Size	Recommend waist range	Belt width 17cm	Belt width 21.5cm
S	75 - 90cm	KLB101	KLB201
M	90 - 105cm	KLB102	KLB202
L	105 - 120cm	KLB103	KLB203
XL	120 - 135cm	KLB104	KLB204
XXL	135 - 150cm	KLB105	KLB205

Not sure how to measure for a belt?
Watch this video



How to use the KomfortLite belt



1

You may find it helpful to start with standing against a wall, with the belt between you and the wall to help hold it in place.



2

Grip the belt on both ends. First, bring the end with the loop pad (the fuzzy pad) around your belly, and place the end with the hook pad (the scratchy pad) over the top.



3

Adjust until you feel comfortable with where the belt is positioned.

TIP:

Ensure the belt is not too tight by checking you can still fit your hand between the belt and your pouch

If you find it hard to grip the belt, use the **new pockets** to help you!



1

Place your hand in the outside hand pocket, hold it against your body to keep it in place.



1B

TIP: It may be helpful to stand against a wall, with the belt between you and the wall to help hold it in place.



2

Place your other hand in the inside hand pocket, this pocket is open on both ends for left or right hands.



3

Bring your hand around your belly, and position it to a position to where you feel most comfortable and snug around your core.



4


Place the first hand over the second hand and adjust as needed.

Tips and resources

View now




Explore our playlist of helpful videos with steps to follow and advice for measuring and applying your support garments.



For active work, sport, or a developed hernia, try the **KoolKnit Support Belt**.




Unsure if this is the right product? **Speak to your STN** about how Support Garments may help you.



View now

For the perfect fit, get measured by your STN. No STN access? Watch the **'How to measure'** video



Worried about Parastomal Hernias?

Ask us for your free copy of our Support Garment Guide which includes information about parastomal hernia risk reduction or if you already have a hernia, how to support and manage it with a support garment.



Just had surgery?

General Guidelines available for returning to normal activities after surgery - written by an Exercise Physiologist!

To find these resources online or to view the full range of support garments, visit omnigon.com.au



OMNIGON

AU 1800 819 274 info@omnigon.com.au
NZ 0800 440 027 www.omnigon.com.au